



ResMed

ResMed Confidential Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way.
© ResMed Ltd 2016.

DOC NO	1031823
PAGE	1 of 2
REV	1

SPECIFICATION – PRINTED MATERIALS

Rev	Change Note	Date document drafted	Document prepared by (Your name)	Document Checked by (Person who will check this document)
1	mC0515	15/9/16	Daniel O'Beirne (or delegate)	Mark Ryan (or delegate)

AirFit F20 Fitting Temp Web mm Eng



ResMed

Actual Size: 140 mm (Width)
x 175 mm (Height)

Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirFit™ F20


Nasal mask
Nasal mask **for Her**

Rx Only  0123



AirFit™ F20

Full face mask
Full face mask **for Her**



© 2016 ResMed Ltd. All Rights Reserved. Specifications may change without notice.
Product availability may vary across regions.

Rx Only  0123



ResMed

Sample size check

10 mm 20 mm 30 mm 40 mm 50 mm 60 mm 70 mm 80 mm 90 mm 100 mm 110 mm 120 mm 130 mm 140 mm 150 mm 160 mm 170 mm